


术前

术后

6个月随访

12个月随访

覆膜支架规格说明图




产品规格范围表

产品规格	规格范围
主体近端直径D1(mm)	26 / 28 / 30 / 32 / 34 / 36 / 38 / 40 / 42 / 44
主体远端直径D2(mm)	20 / 22 / 24 / 26 / 28 / 30 / 32 / 34 / 36 / 38 / 40 / 42 / 44
侧支远端直径D3(mm)	6 / 8 / 10 / 12 / 14
外管外径(Fr)	24
主体长度L1(mm)	60 / 65 / 70 / 75 / 80 / 90 / 100 / 110 / 120 / 130 / 140 / 150 / 160 / 170 / 180 / 190 / 200 / 210
侧支长度L2(mm)	25 / 30 / 35 / 40 / 45
侧支后移长度L3(mm)	5 / 10 / 15 / 20 / 25 / 30



地址: 上海市·浦东新区·康新公路3399弄·1号楼 (上海国际医学园区·医创园)  
 网址: www.endovastec.com  
 电话: (86) (21) 38139300 传真: (86) (21) 33750026 邮编: 201318





- 累及LSA的降主动脉夹层
- 需要将覆膜支架锚定在 ZONE2和ZONE3区域的降主动脉夹层

# 弓部一体，重建分支



# 双入路导入设计，定位精准可靠

多维度摆动角度  
适应不同的分支走向




分支一体化  
无内漏和移位的风险


锥度设计  
适应主动脉夹层的解剖形态



5-30mm分支后移长度




分支根部定位环  
确保侧支通畅圆整




电抛光镍钛丝  
增强耐腐蚀性



高致密覆膜材料  
降低IV型内漏




近端多重覆膜小波段  
提升密封效果




分支入路


主体入路




亲水涂层




软鞘跨弓  
刚柔并济




侧支牵拉  
精确定位



拉线释放  
快速精准



释放分支鞘  
防拉脱设计



弓部一体  
重建分支